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**1. Rights and responsibilities**

It is the responsibility of the customer to carry out the Electrical (including the appropriate earthing), Multimedia, Pneumatic, Water supply and Waste installations, to the connection point of the KaVo product, using the services of a duly qualified person, according to the installation plans supplied by KaVo. The maintenance of the installation cited in the previous item is the responsibility of the owner and/or user of the product. It is not the responsibility of the authorised KaVo technician to measure or check the quality of the earthing (nor the rest of the installation). If the technician finds a fault in the earth wiring, the product will not be installed. The purchaser must pay for the travelling and overstay expenses of the technical staff involved with the installation service call outside the urban perimeter of the city where the authorised service is located. As soon as these instructions are received, the customer should contact the authorised service for clarification.

**Ground characteristics**

The quality of the ground structure must match the load capacity buildings according to the current standards of the site. The floor must be leveled. The chair does not require fastening to the floor. We recommend fastening the dental chair to the floor in the cases where the floor has the following features: wood, hollow blocks, depression, protuberance and unevenness.

**ATTENTION!**

For cases of preparation of a location where there was previously other equipment, it must be re-evaluated according to the standards requested in this document.

**2. Information on the installation for the purchaser****2.1 Electrical Installation**


The distribution box has a connector where the power supply wires are connected. The purchaser must provide the electrical installation in accordance with this instructions. The electrical installation prior to the installation of the product must be done in accordance with the standards in force at the site on "Low voltage electrical installations". The wires of the power supply must have their capacity calculated taking into consideration the consumption of the product, the length of the wires and the maximum voltage drop. Wiring of the trip switch is until the product: 2.5 mm<sup>2</sup>.

Maximum Consumption: 412VA


All KaVo products must be earthed. The quality of the earthing implies the safety of the dentist, his or her auxiliary and the patient. It is recommended that the supply cables are protected by a trip switch 5A. The earthing (as well as the rest of the electrical installation which supplies the product) must be inspected, measured and corrected periodically, as it is subject to degeneration over time. The earth wire must be identified with the green yellow colouring. The neutral wire must be identified with the light blue colouring. It is recommended that the supply cables are protected by a trip switch.

**Monitor Support kit**


The cable signal with the video HDMI (Computer) must have female connector, 60cm out of electrical conduit. The other end (Computer), must have the output of the video signal in perfect condition. The electrical conduit transmission should not have knees 90°, but curves 90°. The KaVo assumes no liability for damage from poor quality of the installation.

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

## English

## 2.7 Area necessary for installation. Positioning of the connection point between equipment and the supply (Model version right-handed)

The measures X and Y are guides and can be altered by the dentist so as to meet the needs.

Guidance values for a patient of 1.76 m in height, with the seat of the chair at a height approximate 0.58 m at its rear part and the backrest in the lying position.

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## 2.2 Pneumatic installation (air)

Compliance with the values below indicated is an indispensable requirement for the perfect functioning of the KaVo equipment. The air supplied to the equipment must be dry, and free of impurities and oil. It must be supplied at a pressure between 5,0 and 5,5 bar. The outflow must be at least 100 liters per minute at 5.5 bar. The air supply branch must have an air filter with a mesh not greater than 25µm, as per the standard ISO 7494-2. The shape and positioning of the tubing must be in accordance with the instructions contained on pages 4, 5, 6 and 7. In order to execute maintenance if there are two or more items of equipment linked to the supply branch, the buyer will be able to adapt an air register in the tubing which is waiting for maintenance. Before installing equipment, see the local legislation concerning air supply, mainly as regards its quality. The main supply network must have a regulator at the input that avoids over-pressure in the components and system.

## 2.3 Water installation

Compliance with the values below indicated is an indispensable requirement for the perfect functioning of the KaVo equipment. The installation and positioning of the tubing must be in accordance with the instructions given in pages 4, 5, 6 and 7. For the purposes of maintenance if there are two or more sets of equipment connected to the same supply branch, the purchaser may adapt one water tap in the wait tubing.

As a safety measure the supply branch must have a general tap which can interrupt the supply when the equipment is not in use.

Before installing the equipment contact the water supply company and find out about the local dispositions concerning the water supply.

It must be free of impurities. If required, pre-connect a water filter.

Quality of the water: drinkable water  
Hardness of the water: 1.50–2.14 mmol/l  
~8.4–12 dH  
Value –pH: 7.2–7.8  
Water filter: 80 m  
Supply pressure: 2.5–6.0 bar  
Supply volume: 3.5 l/min

If exceeding the maximum value, it will increase the tendency of calcareous deposits. If the value is below the minimum, algae is formed.

The main supply network must have a regulator at the input that avoids over-pressure in the components and system.

## 2.4 Waste

The form and positioning of the tubing must be in accordance with the instructions given in pages 4, 5, 6 and 7. The minimum drop is 10mm/m. Observe the standards in force for the locality on “Outflow of used water”.

Nominal diameter of the piping: 40mm

## 2.5 Suction

The form and positioning of the tubing must be in accordance with the instructions given in ages 4, 5, 6 and 7.

Nominal diameter of the piping: 40mm.

In order to avoid deficiencies in the suction, do not use sleeves in the piping or any bends.

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2.6 Pattern for positioning the tubes at the connection point between the equipment and the supply.  
Tubes under the floor.

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

Only to cases of the fixing the floor

A - Equipment fixing points of to the floor. Do not pass tubes under this area.

Referência para posicionamento do ponto de ligação entre equipamento e rede  
Referencia para la posición del punto de conexiones entre el equipo y la red  
Reference for positioning of connection point between equipment and the supply